

## **Skedco HMH SKED® Strap Kit Patient Care | Flexible Litters and Accessories**



The HMH Sked® is designed for use in mass casualty incidents or Haz-Mat incidents. Patients can be rolled or dragged onto the HMH Sked® rather than lifted. This creates less rescuer fatigue and requires fewer rescuers to be exposed in the Hot Zone at a Haz-Mat incident. The HMH Sked® can be easily dragged by personnel wearing level A protection or even pulled by a rope connected outside the Hot Zone.

Key Features: No grommets to trap contaminants • Stores flat in stacks of up to 200 • Ideal for use in mass casualty incidents • Polypropylene straps are chemical resistant • 2" Side release buckles • Not for lifting patients

**Additional Information** 

Shipping weight | 2 lbs Brand | <u>Skedco</u>

**Product Code** 

SK252

